


## Product Datasheet

### PLAT antibody (orb167487)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PLAT                             |
| <b>Reactivity</b>          | Human                                                          |
| <b>Conjugation</b>         | Unconjugated                                                   |
| <b>Tested Applications</b> | WB                                                             |
| <b>Immunogen</b>           | Recombinant fusion protein of human tPA/tPA/PLAT(NP_000921.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.             |
| <b>Concentration</b>       | 1 mg/ml                                                        |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                    |
| <b>Note</b>                | For research use only                                          |
| <b>Isotype</b>             | IgG                                                            |
| <b>Clonality</b>           | Polyclonal                                                     |
| <b>MW</b>                  | 63kDa                                                          |
| <b>Uniprot ID</b>          | <b>P00750</b>                                                  |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                |



Western blot analysis of extracts of var...