



Product Datasheet

TET2 antibody (orb167514)


Description	Rabbit polyclonal antibody to TET2.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human TET2(NP_001120680.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	280KDa
Uniprot ID	Q6N021
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A431...



Immunohistochemical staining of rat lung...



Immunohistochemical staining of human lu...