


Product Datasheet

EGFR (phospho-Y1092) antibody (orb167531)

Description	Rabbit polyclonal antibody to EGFR
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around Y1068 of human EGFR(NP_005219.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	185kDa
Uniprot ID	P00533
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...