

Product Datasheet

Human Cadherin-17 / CDH17 Protein, Mouse IgG1 Fc Tag (orb1675334)

Description

Human Cadherin-17 Protein,

from human 293 cells (HEK293). It contains...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Cadherin-17

Preservatives

20 mM Tris, 150 mM NaCl, pH8.0

Form/Appearance

Powder

Storage

-20°C

Tag

C-mFc(mIgG1)

Note

For research use only

Application notes

This protein carries a mouse IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 111.2 kDa. The protein migrates as 130-135 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Gln 23 - Met 787

Purity

90%

Source

Human Cadherin-17 Protein, Mouse IgG1 Fc Tag (orb1675334) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Met 787 (Accession # Q12864-1).

MW

111.2 kDa

Expiration Date

6 months from date of receipt.