

Product Datasheet

Mouse Apolipoprotein E Protein, His Tag (MALS & SPR verified) (orb1675335)

Description

APOE is an apolipoprotein, a

particles, that mainly functions in lipo...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Apolipoprotein E

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 35.9 kDa. The protein migrates as 14 kDa and 35 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence

Glu 19 - Gln 311

Purity

90%

Source

Mouse Apolipoprotein E Protein, His Tag (orb1675335) is expressed from human 293 cells (HEK293). It contains AA Glu 19 - Gln 311 (Accession # P08226).

MW

35.9 kDa

Expiration Date

6 months from date of receipt.