

Product Datasheet

Human GDNF / ATF / hGDNF Protein, premium grade (orb1675340)

Description Glial cell line-derived neurotrophic factor (GDNF) is a neurotrophic factor that enhances survival a...

Reactivity Human

Endotoxins 0.1 EU per µg

Conjugation Unconjugated

Target GDNF

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag Native

Note For research use only


Application notes This protein carries no "tag". The protein has a calculated MW of 15.1 kDa. The protein migrates as 15-22 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity 90%

Source Human GDNF Protein, premium grade (orb1675340) is expressed from human 293 cells (HEK293). It contains AA Ser 78 - Ile 211 (Accession # P39905-1).

MW 15.1 kDa

Expiration Date 6 months from date of receipt.



Immobilized
Human
GDNF
Protein,
premium ...