

Product Datasheet

**Human BDNF / Abrineurin Protein, premium grade
(orb1675341)**

Description

Human BDNF, premium grade (BDF-H5219) is

contains AA His...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

BDNF

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;
-70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 13.5 kDa. The protein migrates as 15-17 kDa and 19 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human BDNF, premium grade (orb1675341) is expressed from human 293 cells (HEK293). It contains AA His 129 - Arg 247 (Accession # P23560-1).

MW


13.5 kDa

Uniprot ID

P23560

Expiration Date

6 months from date of receipt.



Immobilized
Human
BDNF,
premium
grade (C...