

Product Datasheet

Human KIR2DL5 / CD158f Protein, Fc Tag (orb1675345)

Description

Killer cell immunoglobulin-like receptor

immunoglobulin superfam...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

KIR2DL5A

Preservatives

50 mM Tris,100 mM Glycine,25 mM Arg,150 mM NaCl,pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50 kDa. The protein migrates as 63-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

His 22 - His 240

Purity

95%

Source

Human KIR2DL5 / CD158f Protein, Fc Tag (orb1675345) is expressed from human 293 cells (HEK293). It contains AA His 22 - His 240 (Accession # Q8N109-1).

MW

50 kDa

Expiration Date

6 months from date of receipt.