


Product Datasheet

Anti-SARS-CoV-2 (2019-nCoV) Nucleocapsid (orb1675357)

Description	Anti-SARS-CoV-2 (2019-nCoV) Nucleocapsid
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419 aa)
Target	SARS-CoV-2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	46
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-SARS-Co...