


Product Datasheet

Anti-MIST1 (orb1675367)

Description	Anti-MIST1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombined human BHLHA15 protein (1-189/189 aa)
Target	BHLHA15
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	21
Uniprot ID	Q7RTS1
Dilution Range	WB=1:100-1000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human K562 cell; Lane...