


Product Datasheet

Anti-Interleukin-1 receptor-like 1/ST2 (orb1675368)

Description	Anti-Interleukin-1 receptor-like 1/ST2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human human Interleukin-1 receptor-like 1/ST2 protein (19-328/369 aa)
Target	IL1RL1
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	61
Uniprot ID	Q01638
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Recombinant human ST2 protein