


Product Datasheet

Anti-Alpha 1 microglobulin (orb1675377)

Description	Anti-Alpha 1 microglobulin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Alpha 1 Microglobulin (20-203/352 aa)
Target	AMBP
Preservatives	PBS
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	37
Uniprot ID	P02760
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Recombinant human Al...