

Product Datasheet

Mouse Complement C5a ELISA Kit (orb1675403)

Description	Mouse Complement C5a ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: C5. Gene_full_name: Complement C5. Protein_name: Complement C5. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: epigenetics and nuclear signaling transcription domain families hlh leucine zipper hlh immunology innate immunity mast cells chip"ing antibodies. Sequence_similarities: Contains 1 anaphylatoxin-like domain.
Expiration Date	Please enquire.