

Product Datasheet

Mouse Cd69/Early activation antigen CD69 ELISA Kit (orb1675603)

Description	Mouse Cd69/Early activation antigen CD69 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78-5000 pg/mL. Sensitivity: less than 39.4 pg/mL. Assay Type: Sandwich. Assay Time: 3-5 h. Detection wavelength: 450 nm. Gene_full_name: Cd69. Sample Volume: 100 µL. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Expiration Date	Please enquire.