

## Product Datasheet

### Fish Thyroxine,T4 ELISA Kit (orb1675731)

|                            |                                                                                                                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Fish Thyroxine,T4 ELISA Kit                                                                                                                                                                  |
| <b>Reactivity</b>          | Fish                                                                                                                                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                        |
| <b>Note</b>                | For research use only                                                                                                                                                                        |
| <b>Application notes</b>   | Detect Range: 20 ng/mL-320 ng/mL.<br>Sensitivity: 20 ng/mL. Assay Time: 1-5 h.<br>Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma and tissue homogenates. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                              |