

Product Datasheet

Chicken Testosterone (T) ELISA Kit (orb1675734)

Description	Chicken Testosterone (T) ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0 ng/mL-20 ng/mL. Sensitivity: 0.024 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.