

Product Datasheet

Cleaved PARP (Asp214) Rabbit Antibody (orb1675779)

Description	Cleaved PARP (Asp214) Rabbit Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Cleaved PARP (Asp214)(Asp214).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89KDa
Uniprot ID	P09874
Dilution Range	WB,1:1000 - 1:2000
Expiration Date	12 months from date of receipt.