

Product Datasheet

ACE2 Rabbit Antibody (orb1675812)

Description	ACE2 Rabbit Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human ACE2(NP_068576.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63kDa/92kDa
Uniprot ID	Q9BYF1
Dilution Range	IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.