

Product Datasheet

Catalase Rabbit Antibody (orb1675822)

Description	Catalase Rabbit Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Recombinant fusion protein of human Catalase(NP_001743.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59kDa
Uniprot ID	P04040
Dilution Range	IP,50 - 1:200
Expiration Date	12 months from date of receipt.