

Product Datasheet

DNM1L Rabbit Antibody (orb1675871)

Description	DNM1L Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human DNM1L(NP_036192.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa/78kDa/79kDa/80kDa/81kDa/82kDa/83kDa
Uniprot ID	000429
Dilution Range	IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.