

Product Datasheet

LC3B Rabbit Antibody (orb1675891)

Description	LC3B Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human LC3B(NP_073729.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	Q9GZQ8
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:100 IF/ICC,1:50 - 1:100 IP,1:20 - 1:50
Expiration Date	12 months from date of receipt.