

Product Datasheet

VDAC1 / Porin Rabbit Antibody (orb1675897)

Description	VDAC1 / Porin Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human VDAC1 / Porin(NP_003365.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32KDa
Uniprot ID	P21796
Dilution Range	WB,1:500 - 1:2000 IHC,1:100 - 1:200
Expiration Date	12 months from date of receipt.