

Product Datasheet

Phospho-p38 MAPK-T180/Y182 Rabbit Antibody (orb1675905)

Description	Phospho-p38 MAPK-T180/Y182 Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around T180 & Y182 of human p38 MAPK(NP_620581.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	35KDa
Uniprot ID	Q16539 , O15264 , Q15759 , P53778
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.