

Product Datasheet

Kallistatin (SERPINA4) Rabbit Antibody (orb1675921)

Description	Kallistatin (SERPINA4) Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Kallistatin (SERPINA4)(NP_006206.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49KDa
Uniprot ID	P29622
Dilution Range	WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:100
Expiration Date	12 months from date of receipt.