

Product Datasheet

IRAK4 Rabbit Antibody (orb1676028)

Description	IRAK4 Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human IRAK4(NP_057207.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37kDa/51kDa
Uniprot ID	Q9NWZ3
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.