

Product Datasheet

M6PR Rabbit Antibody (orb1676103)

Description	M6PR Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human M6PR(NP_002346.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	P20645
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.