

Product Datasheet

δ -Catenin/p120 Catenin Rabbit Antibody (orb1676129)

Description	δ -Catenin/p120 Catenin Rabbit Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human δ -Catenin/p120 Catenin(NP_001078938.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65-72kDa/90-108kDa
Uniprot ID	O60716
Dilution Range	WB,1:500 - 1:2000 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.