

Product Datasheet

VPS18 Rabbit Antibody (orb1676183)

Description	VPS18 Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human VPS18VPS18
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110kDa
Uniprot ID	Q9P253
Dilution Range	WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:200
Expiration Date	12 months from date of receipt.