

Product Datasheet

P4HB Rabbit Antibody (orb1676199)

Description	P4HB Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human P4HB(NP_000909.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P07237
Dilution Range	WB,1:500 - 1:1000 IHC,1:50 - 1:200 IF/ICC,1:50 - 1:100
Expiration Date	12 months from date of receipt.