

Product Datasheet

Phospho-PKC mu-S910 Rabbit Antibody (orb1676203)

Description	Phospho-PKC mu-S910 Rabbit Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S910 of human PKC mu(NP_002733.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	101kDa
Uniprot ID	Q15139
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.