


Product Datasheet

Human B7-H2/ICOSL ELISA Kit (orb1676355)

Description	Human B7-H2/ICOSL ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D19-S258
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	075144
Sensitivity	300 pg/ml
Expiration Date	Please enquire.



Human
ICOSLG
PicoKine
ELISA Kit
standard...