

Product Datasheet

ARS2/SRRT Antibody (orb1676512)

Description

ARS2/SRRT Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human ARS2/SRRT recombinant protein (Position: H125-Q793).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Clonality

Polyclonal

MW


125 kDa

Uniprot ID
Q9BXP5
Dilution Range


Western blot, 0.25-0.5 µg/ml, Human Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5 µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of ARS2/SRRT using anti-ARS...



IHC analysis of ARS2/SRRT using anti-ARS...



IHC analysis of ARS2/SRRT using anti-ARS...