


Product Datasheet

Mouse Nepmucin/CD300g ELISA Kit (orb1676581)

Description	Mouse Nepmucin/CD300g ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L19-I284
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q1ERP8
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Nepmucin
PicoKine
ELISA Kit
standa...