

Product Datasheet

EV71 Antibody (orb1676673)

Description	Enterovirus EV71 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, EIA, LF
Target	EV71
Preservatives	0.05% Proclin 300
Form/Appearance	Supplied as a 0.22 µm filtered solution of PBS, pH 7.4
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Application notes	Use: Conjugate
Isotype	IgG
Clonality	Recombinant
Purity	≥90% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.