

Product Datasheet

TSH Antibody (orb1676684)

Description	Thyroid stimulating hormone(TSH) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, EIA, LF
Target	TSH
Preservatives	0.05% Proclin 300
Form/Appearance	Supplied as a 0.22 µm filtered solution of PBS, pH 7.4
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	≥90% as determined by SDS-PAGE
Expiration Date	12 months from date of receipt.