


Product Datasheet

SLC26A6 antibody (orb167675)

Description	Goat polyclonal antibody to SLC26A6.
Species/Host	Goat
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SLC26A6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in Human and Mouse Pancreas lysates (calculated MW of 80.8kDa according to Human NP_602298.2 and 80.5kDa according to Mouse NP_599252.2). Recommended concentration: 0.1-0.3µg/ml. Primary incubation 1 hour at room temperature.
Protein Sequence	RDVAEYSEAKEVR
Clonality	Polyclonal
MW	83.0; 82.8; 80.8; 80.9
Entrez	65010
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Pancreas ...