


Product Datasheet

POU2F2 antibody (orb167679)

Description	Goat polyclonal antibody to POU2F2
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POU2F2 / OCT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Tonsil and Mouse Thymus lysates (calculated MW of 49.5kDa according to NP_002689.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Pelis et al, Am J Physiol Renal
Protein Sequence	KQGLDSPSEHTDTER
Clonality	Polyclonal
MW	49.5
Entrez	5452
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Tonsil ly...