


Product Datasheet

Cntn4 antibody (orb167682)

Description	Goat polyclonal antibody to Cntn4
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Contactin 4 / Big-2 (mouse aa160-172)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in Mouse Brain lysates and not in the same lysates from the KO mice (calculated MW of 113kDa according to NP_001103219.1). Recommended concentration: 0.5-1.5µg/ml.
Protein Sequence	RRADGKPIARK
Clonality	Polyclonal
MW	113; 78.4
Entrez	152330
Dilution Range	ELISA: 1:128000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Olfactory...