


Product Datasheet

DDIT3 antibody (orb167683)

Description	Goat polyclonal antibody to DDIT3
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CCAAT/enhancer binding protein zeta (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 140kDa band observed in Mouse fetal Brain lysates (calculated MW of 120kDa according to NP_001019977.1). Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	DNASFKQLKWEAERD
Clonality	Polyclonal
MW	120
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse fetal Brain...