


Product Datasheet

PYROXD1 antibody (orb167686)

Description	Goat polyclonal antibody to PYROXD1.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	PYROXD1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Small Intestine. Recommended concentration: 5µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in U2OS cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of U2OS cells. Recommended concentration: 10ug/ml.
Protein Sequence	TCAEQLATHFPSED
Clonality	Polyclonal
MW	55.8
Entrez	79912
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Colorectal cancer tissue.