



Product Datasheet

LIF antibody (orb167688)

Description	Goat polyclonal antibody to LIF
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	LIF (aa28-39)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 23kDa band observed in Mouse and Rat Colon lysates (calculated MW of 22.3kDa according to Mouse NP_032527.1 and 22.1kDa according to Rat NP_071532.2). Recommended concentration: 0.5-1.5µg/ml. Primary incubation 1 hour at room temperature.
Protein Sequence	TPVNATCAIRHP
Clonality	Polyclonal
MW	22
Entrez	3976
Dilution Range	ELISA: 1:64000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Colon lys...



Western blot analysis of Mouse(lane 1) a...