



Product Datasheet

Cd97 antibody (orb167690)

Description	Goat polyclonal antibody to Cd97
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Cd97 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in Mouse Thymus and Spleen lysates (calculated MW of 80.2kDa according to NP_001156503.1). The beta chain is detected as a ~30kDa band in Mouse Colon (transgenic overexpressing Cd76 in intestines) and not in the knock-out. Reco
Protein Sequence	TSQTRALRSSES
Clonality	Polyclonal
MW	90.4; 80.4; 85.6; 80.2
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Thymus ly...



Western blot analysis of Mouse Colon (wt...