

Product Datasheet

NANOG monoclonal antibody (orb1676944)

Description	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human NANOG. The exact sequence is proprietary.
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 35 kDa
Uniprot ID	Q9H9S0
Dilution Range	WB (1/1000 - 1/2000), IHC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.