

Product Datasheet

PRDM12 monoclonal antibody (orb1676946)

Description	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PRDM12. The exact sequence is proprietary.
Preservatives	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 40 kDa
Uniprot ID	Q9H4Q4
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.