

Product Datasheet

THBS1 antibody (orb167695)

Description

Goat polyclonal antibody to THBS1.

Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Target

Thrombospondin

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: Approx 150kDa band observed in Mouse Kidney lysates (calculated MW of 129kDa according to NP_035710.2). Recommended concentration: 0.3-1µg/ml. Primary incubation 1 hour at room temperature.. IHC: Paraffin embedded Human Spleen. Recommended concentration: 5µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of HeLa and HepG2 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration

Protein Sequence

NRIPESGGDNSVFD

Clonality

Polyclonal

MW

129

Entrez


7057

Dilution Range


ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of HEK293 lysate (...)



Immunohistochemical staining of Human Sp...



Immunohistochemical staining of Human Pl...