

Product Datasheet

p63 monoclonal antibody (orb1676951)

Description	Mouse IgM. Supplied in crude ascites with 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human p63. The exact sequence is proprietary.
Preservatives	Mouse IgM. Supplied in crude ascites with 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 70 kDa
Uniprot ID	Q9H3D4
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.