


Product Datasheet

SRCAP antibody (orb167696)

Description	Goat polyclonal antibody to SRCAP
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SRCAP (aa309-322)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Skin shows nuclear staining in some keratinocytes. Recommended concentration: 5-10µg/ml.
Protein Sequence	QEGNDAEAQRREIE
Clonality	Polyclonal
MW	344
Entrez	10847
Dilution Range	ELISA: 1:128000, IHC-P: 5-10 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of Human
Sk...