

Product Datasheet

CRMP-5 monoclonal antibody (orb1676967)

Description	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human CRMP-5. The exact sequence is proprietary.
Preservatives	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 61 kDa
Uniprot ID	Q9BPU6
Dilution Range	WB (1/1000 - 1/2000), IF/ICC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.