


Product Datasheet

CHUK antibody (orb167700)

Description	Goat polyclonal antibody to CHUK.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CHUK
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 84.6kDa according to NP_001269.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	EHDHSLSCVVPQD
Clonality	Polyclonal
MW	84.6
Entrez	1147
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa nuclear lysates.