



Product Datasheet

APOD antibody (orb167709)


Description	Goat polyclonal antibody to APOD.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	APOD
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 25kDa band observed in Human Liver and Kidney lysates (calculated MW of 21.3kDa according to NP_001638.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.3-1µg/ml.. IHC: Paraffin embedded Human Kidney and Adrenal Gland. Recommended concentration: 5µg/ml.
Protein Sequence	EIEKIPTTFENGRC
Clonality	Polyclonal
MW	21,3
Entrez	347
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Ki...



Immunohistochemical staining of Human Ad...